ARMY INST OF DENTAL RESEARCH WASHINGTON DC
PRE-CEMENTATION STABILIZING OF ETCH-ON FIXED PROSTHESES, (U)
WAR 82 L LORTON, D DENUCCI
UNCLASSIFIED


M.

Harl MARK





END PATE FILMED 5 -82



MICROCOPY RESOLUTION TEST CHART

A112993

REPORT DOCUMENTATION PAGE		READ INSTRUCTIONS BEFORE COMPLETING FORM
REPORT NUMBER	2. GOVT ACCESSION NO.	3. RECIPIENT'S CATALOG NUMBER
TITLE (and Subtitie)		5. TYPE OF REPORT & PERIOD COVERED
Pre-Cementation Stabilizing of Etch-On Fixed Prostheses		N/A
		6. PERFORMING ORG. REPORT NUMBER
AUTHOR(e)		8. CONTRACT OR GRANT NUMBER(a)
Lewis Lorton and D. DeNucci		N/A
PERFORMING ORGANIZATION NAME AND ADDRESS		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS
JS Army Institute of Dental Research Valter Reed Army Medical Center		N/A
CONTROLLING OFFICE NAME AND ADDRESS		12. REPORT DATE
		25 March 1982
		13 NUMBER OF PAGES
MONITORING AGENCY NAME & ADDRESS(If different from Controlling Office)		15. SECURITY CLASS (at this report)
J.S. Army Medical Research & Deve	elopment Command	
HQDA-IS		UNCLASSIFIED
Fort Detrick, MD 21701		15a. DECLASSIFICATION DOWNGRADING SCHEDULE


This document has been approved for public release and sale; its distribution is unlimited.

- 17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, If different from Report)
- 18. SUPPLEMENTARY NOTES
- 19. KEY WORDS (Continue on reverse side if necessary and identify by block number)

Fixed Prostheses

20. ABSTRACT (Continue on reverse elde if necessary and identity by block number)

N/A



EDITION OF 1 NOV 65 IS OBSOLETE

UNCLASSIFIED

Pre-cementation stabilization of etch-on fixed prostheses

*L. Lorton, D.D.S., M.S.D. **D. DeNucci, D.D.S., M.S.

The etch-on fixed prostheses is particularly useful in the replacement of missing lower anterior teeth where there is also some need for periodontal splinting.

There is an obstacle to accurate cementation that is inherent in the design of these "minimum-framework" prostheses. Although the framework may have a positive seat when there is no intervening composite, as soon as the cementing medium is introduced, all tactile guides to seat the framework are lost.

Temporary incisal rests. cast with the framework, will provide a positive orientation for the seated framework. These incisal rests can often be placed so as not to interfere with occlusal adjustments of the uncemented restoration. To facilitate occlusal adjustment or their subsequent removal, the rests may be ground to a minimum thickness.

After the composite has set, the rests may be easily ground off, leaving the correctly oriented prosthesis.

The opinions contained herein are the private views of the authors and are not to be construed as official or as reflecting the views of the Dept. of the Army or the Dept. of the Defense

Reprint requests to: LTC Lewis Lorton Division of Clinical Operations US Army Institute of Dental Research Walter Reed Army Medical Center Washington, DC 20012

^{*}LTC (DC) USA, Assistant Dental Research Coordinator, US Army Institute of Dental Research, Walter Reed Army Medical Center, Washington, DC 20012 **LTC (DC) USA, Senior Periodontic Resident, US Army Institute of Dental Research, Walter Reed Army Medical Center, Washington, DC 20012

FIGURE 1 Framework in place showing incisal rests.

FIGURE 2 Incisal rests at midpoint of onlayed tooth.



